





# Anti-Decubitus Systems

have been designed to ensure effective prevention and therapy of all stages of pressure ulcers. HICO's range of decubitus treatment products offers your patients soothing pressure and pain relief. The washable and disinfectable alternating pressure mattresses provide optimum patient comfort. Thanks to the breathable, non-slip mattress cover, shear forces are minimized - thus reducing the risk of bedsores.

#### **Properties**

- Suitable for patients with a body weight of 25 to 250 kg (depending on the product)
- High-quality 2- or 3-chamber air mattresses
- Low-noise and -vibration devices, cycle times adjustable between 6 and 30 min. (depending on the product)
- Transport mode and static pressure for soft lying comfort
- Automatic pressure control and fault diagnosis
- CPR valve for quick emergency air release (depending on the product)

#### 25 - 250 kg

#### **HICO-DECUBIMAT 390**

Powerful alternating pressure system for patients with advanced stages of pressure ulcer



#### 40 - 160 kg

#### HICO-DECUBIMAT 370

For the pain-reducing positioning of patients suffering from mediumstage pressure ulcers



#### 40 – 120 kg

#### **HICO-DECUBIMAT 365**

For the prevention of pressure ulcers in nursing homes



#### 20 – 120 kg

#### **HICO-FOAM**

Non-energetic anti-decubitus overlay mattress for providing pressure relief



## Medical-technical expertise for more than 70 years

Founded by Hans Hirtz in 1945, the Cologne-based company Hirtz & Co. KG is run in the fourth generation by the Hirtz family. Today, the company stands for expertise in medical technology sold under the brand name HICO Medical Systems. What started as a wholesale business for medical products and devices has developed - thanks to in-house innovations - into a major medtech company. The product portfolio comprises hypo- and hyperthermia systems, decubitus prevention and treatment solutions as well as inhalation devices

Since 1997, Hirtz & Co. KG has been certified according to ISO standards and therefore guarantees excellent quality. Today, more than 70 years after the family business was founded, brand products of HICO Medical Systems generate more than 50 % of the company's sales outside Germany. To ensure our future in the market with leading medtech products, Hirtz & Co. KG places great value on well-established development partnerships with leading university hospitals.

## **HICO-DECUBIMAT 390**

Powerful alternating pressure system for patients with advanced stages of pressure ulcer



HICO-DECUBIMAT 390 combines the advantages of our most powerful alternating pressure pump with those of the 3-chamber deep-cell mattress. The result: effective pressure relief for patients with a body weight of up to 250 kg.

#### **Properties**

- Periods of "zero pressure" and optimum reduction of contact pressure for patients with a weight of 25 to 250 kg
- 3-chamber deep-cell mattress with CPR valve and static head cell
- Transport mode and static pressure for soft support
- High lying comfort thanks to breathable, non-slip mattress cover (PU-coated, bi-elastic nylon)
- Extremely low noise and vibration levels



#### Device

REF 340001

Dimensions W x H x D approx. 40 x 35 x 15 cm

Noise emission approx. 30 dB(A) (1 m)

#### **Mattress**

RFF 340040

3-chamber system

Dimensions L x W x H approx. 200 x 90 x 21 cm

Weight approx. 9.7 kg



### **HICO-DECUBIMAT 370**



For the pain-reducing positioning of patients suffering from medium-stage pressure ulcers



HICO-DECUBIMAT 370 features a large-cell, 2-chamber alternating pressure mattress for the pain-relieving positioning of patients (up to 160 kg weight) suffering from pressure ulcers.

#### **Properties**

- Patient weight adjustable within a range of 40 to 160 kg
- 2-chamber large-cell mattress with CPR valve and static head cell
- Transport mode and static pressure for soft support
- High lying comfort thanks to breathable, non-slip mattress cover (PU-coated, bi-elastic nylon)
- Low noise and vibration levels



#### **Device**

REF 320003

Dimensions W x H x D approx.  $27 \times 25 \times 12 \text{ cm}$ 

Noise emission approx. 28 dB(A) (1 m)

#### **Mattress**

REF 350020

2-chamber system

Dimensions L x W x H approx. 200 x 90 x 13 cm

Weight approx. 6.2 kg



## **HICO-DECUBIMAT 365**

For the prevention of pressure ulcers in nursing homes



Especially in care facilities, HICO-DECUBIMAT 365 is used for the prevention of pressure ulcers and for the pain-relieving positioning of patients (up to 120 kg weight) suffering from early-stage bedsores.

#### **Properties**

- Good pressure relief for patients with a body weight of 40 to 120 kg
- Alternating pressure support mattress made of skin-friendly vinyl
- Airflow system for patient ventilation through very small holes (laser perforation)



#### **Device**

REF 310003

Dimensions W x H x D approx. 27 x 25 x 12 cm

Noise emission approx. 23 dB(A) (1 m)

#### **Mattress**

REF 310500

2-chamber system

Dimensions L x W x H approx. 200 x 90 x 7 cm

Weight approx. 2 kg



## **HICO-FOAM**



Non-energetic anti-decubitus overlay mattress for providing pressure relief



HICO-FOAM has established itself as a non-energetic, pressure-relieving overlay for patients up to a body weight of 120 kg.

#### **Properties**

- Pressure relief for patients with a body weight of 20 to 120 kg
- 3-layer mattress system made of open-pored MDI cold foam
- Ideal combination of longitudinal and transverse sections with different compression hardnesses to ensure effective pressure relief
- Bi-elastic, water-vapor permeable PU cover for minimizing shear forces



REF 310028

Non-energetic overlay

Dimensions L x W x H approx. 200 x 90 x 15 cm

Weight approx. 11 kg



## **Further product lines**







