

TOOTH PROTECTOR

Indispensable anesthetic aid during intubation



The TOOTH PROTECTOR provides ideal protection of the teeth during intubation. The dental guard fits every size of mouth and cannot be swallowed.

Properties

- Reduces the risk of injury and protects the teeth
- Opens the mouth
- Fits every size of mouth and cannot be swallowed
- Jumps out of the mouth if fitted incorrectly
- Latex-free
- Non-sterile
- For one-time use only - cannot be autoclaved
- Supplied in a carton with 100 single packs



REF 640012

Material
Food-safe TPE



Further product lines



Patient Temperature Management

Precise warming and cooling of your patients during surgery and in intensive care for speeding up recovery



Anti-Decubitus Systems

Soothing pressure relief and extra high patient comfort



Ultrasonic Nebulization

A wide range of different uses to ensure effective inhalation therapy

